

THE LONGTERM COOLING INLAY

Our own design, to achieve a super cooling breakthrough for temperature sensitive medication

NOVACOOOL Temperature Protective Thermobags, are already being used worldwide by many Mega Pharmaceutical Companies, such as Novartis and Sanofi, with great success.

Due to innovative material development, **NOVACOOOL Temperature Protective Thermobags**, can already protect temperature sensitive medication, with just the use of one small 200grm freezer pack, over a constant temperature range of 2 to 8 degree C, for up to 9.5 hours.

This achievement was officially certified for the Austrian Firm NovaLogo, by the Austrian Forschungs Institut in Vienna.

Ongoing research by NovaLogo, to obtain for clients and patients an even longer temperature protection time, without having to change the cooling freezer packs, has now resulted in a phenomenal solution, namely:

THE LONGTERM COOLING INLAY

This newly developed product is a special bag, designed to hold a 200grm bubble wrapped freezer pack and the to be cooled medication. When inserted into a Nova-Logo Thermo Bag, this combination of NOVACOOOL THERMOBAG, coupled with the LONGTERM COOLING INLAY, significantly increases the temperature protection time from 9.5 hours to 17 hours at a constant temperature range of 2-8 degrees C.

Additionally, with the use of a 400grm freezer pack, this combination achieves the unheard of protection time of 28 hours, at a constant temperature range of 2 to 8 degrees C!

With such long lasting protection time, temperature sensitive medication can now be transported, safely cooled, over long distances, such as Europe to Africa, Asia, or North and South America. And all this, without the use of electricity.



NOVACOOOL
THERMOBAG

Step by Step User Information:



- 1 Longterm Cooling Inlay
- 2 NOVACOOOL Thermobag
- 3 Frozen Freezerpack inside open Bubble Wrap Envelope
- 4 Insert FrozenFreezer Pack with opened Bubble Wrap Envelope at bottom of Longterm Inlay and place to be cooled product on top .
- 5 Put LONGTERM COOLING INLAY with contents into open, empty NOVACOOOL Thermobag
- 6 Close Zippers of both the Longterm Cooling Inlay and the NOVACOOOL Thermobag.



Tel. 613-662-3515
info@novalogomedical.com
novalogomedical.com

Address:
7-841 Sydney Street, suite 175
Cornwall, Ontario
Canada, K6H 7L2